

DuPont Qualicon BAX® System

REAL-TIME PCR ASSAY

FOR *CAMPYLOBACTER JEJUNI/COLI/LARI*

PRODUCT DESCRIPTION

Campylobacter is a worldwide public health concern and the leading cause of enteric illness in the US. Ingesting even low doses (less than 500 cells) can cause Campylobacteriosis, with clinical complications that include arthritis and Guillain-Barré syndrome. Although *C. jejuni* is responsible for 90% of *Campylobacter* infections, *C. coli* and *C. lari* are also harmful to humans. These bacteria are carried in the intestinal tracts of animals, especially birds. The foodborne pathogen has been found in unpasteurized milk, meat, poultry, shellfish, fruits and vegetables. Current screening procedures are culture-based, which require a laborious plating method and at least five days to result.



BAX® System Real-Time PCR Assay *Campylobacter jejuni/coli/lari*

Part # D12683449

96 tests per kit

PCR tubes with tablets, optical caps,
protease, lysis buffer

Store at 2-8°C

Stable to expiration date on label

Benefits of the BAX® System Real-Time PCR Assay

- Speed – Same-day results for highly contaminated samples without enrichment, next-day results for enriched samples
- Accuracy – Detects and quantifies all three species in the same sample
- Exceptional sensitivity – Reliably detects 10⁴ cfu/mL
- Ease of use – Tableted reagents reduce operator error
- Closed-cap system avoids amplicon contamination in the lab
- LIMS-compatible electronic data for easy storage, sharing and retrieval

Features

- Less than 90 minutes real-time processing
- 24-48 hours enrichment
- If any species is present, sample is positive; if none is present, sample is negative
- Displays separate quantitative values as cfu/mL for each species in the sample – no interpretation required
- Validated on enriched samples of spiked turkey breast and naturally contaminated chicken carcass rinses
- Sensitivity equivalent to ISO 10272-1:2006(E) culture-based method
- Specificity ≥99%
- Excellent inclusivity/exclusivity

Designed for efficient workflow and reliable results

- 96 tests per kit
- Lysis reagents included
- Can be combined with other BAX® system real-time assays
- Step-by-step directions included in each package

Certifications



AOAC Research Institute
Performance Tested Methodsm
#040702 – Validated on processed
turkey and poultry carcass rinses

NordVal #39

Validations and Approvals

- People's Republic of China AQSIQ
- Russia Rospotrebnadzor



The miracles of science™

Sample preparation



Without enrichment: Rinse poultry carcass in 400 mL buffered peptone water.


With enrichment:

Processed poultry product: Prepare 1:10 dilutions of sample in Bolton broth with supplement.

Carcass rinses: Rinse poultry carcass in 400 mL of buffered peptone water, and add 30 mL rinse to 30 mL of double-strength Bolton broth. Incubate all samples under microaerobic conditions at 42 ±1°C for at least 24 hours. You can extend enrichment to 48 hours if needed.

BAX® system protocol


8:00 Create rack file and warm up cyclor.



8:20 Heat cluster tubes for 20 minutes at 37°C, then 10 minutes at 95°C.




10:30 Review results.




Target	Result	Ct	CFU/ml
<input checked="" type="checkbox"/> C. coli	Positive	32.8	1.3E+06
<input checked="" type="checkbox"/> C. jejuni	Positive	33.7	1.3E+06
<input checked="" type="checkbox"/> C. lari	Positive	29.8	1.4E+06
<input checked="" type="checkbox"/> IPC	Positive	34.4	


8:05 Mix protease with lysis buffer and transfer 200 µL of lysis reagent to cluster tubes.



8:50 Cool cluster tubes for 5 minutes in cooling block, then transfer 30 µL to PCR tubes in cooling block.



8:10 Transfer 5-µL samples to cluster tubes.



9:00 Place sealed PCR tubes in cyclor and run program.



DuPont Qualicon
 ESL Bldg 400, PO Box 80400
 Wilmington, DE 19880-0400 USA
 Tel: 800-863-6842 or 302-695-5300
 Fax: 302-695-5301
 Europe: 00 800 3876 6838
 Singapore Tel: +65 6586 3635
 www.qualicon.com



The miracles of science™